


# Robust Plus

Aurelia® Robust Plus™ Extended  
Cuff Blue Nitrile Powder Free  
Examination Gloves



**Extended Cuff**  
**Maximum Protection**

- 100% latex free, blue color
- 5 mil thickness for improved puncture resistance
- Micro-textured finish for improved gripping action
- 12" Extended Cuff




**AURELIA**  
Sign of Comfort™

[www.aureliagloves.com](http://www.aureliagloves.com) | 1-877-AURELIA (1-877-287-3542)

# Robust Plus

## Product Description

Aurelia® Robust Plus™ Extended Cuff Blue Nitrile (Non Latex) Powder Free Examination Gloves are made from 100% Nitrile (Non Latex), ambidextrous, blue color, Finger Textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit [www.aureliagloves.com](http://www.aureliagloves.com)



## Aurelia® Robust Plus™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Small	26.1	550
Medium	26.5	550
Large	26.1	595
Extra Large	26.2	580
ASTM D6319-00	14 MPa minimum	500% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Small	24.0	550
Medium	24.0	555
Large	24.0	585
Extra Large	23.7	570
ASTM D6319-00	14 MPa minimum	500% minimum

## DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Small</b> Aurelia Robust Plus ASTM D6319-00	12 / 305 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Medium</b> Aurelia Robust Plus ASTM D6319-00	12 / 305 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Large</b> Aurelia Robust Plus ASTM D6319-00	12 / 305 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
<b>Extra Large</b> Aurelia Robust Plus ASTM D6319-00	12 / 305 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

### ITEM NUMBER

Small	Medium	Large	Extra Large	XX-Large
63886	63887	63888	63889	63880

## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor: