


Dimensions


Aurelia® Dimensions™ Fitted Powder Free
Latex Examination Gloves



**Hand-Specific Comfort
Superb Ergonomics**



- Left and right hand precision-fit gloves for most comfortable fit
- Reduces risk of carpal tunnel syndrome
- Fingertip-textured finish for precise gripping
- Polymer-coated helps gloves slide on with ease, keeping hand smooth



AURELIA®
Sign of Comfort™


www.aureliagloves.com | 1-877-AURELIA (1-877-287-3542)

Dimensions

Product Description

Aurelia® Dimensions™ Fitted

Powder Free Latex Examination Gloves are made from 100% natural Malaysian rubber, Hand-Specific, creamy white color, Fingertip textured and non-sterile. Aurelia gloves conform and are tested to the highest international standards. For more information visit www.aureliagloves.com



Aurelia® Dimensions™ Specifications

TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	970
Small	31.1	970
Medium	31.0	970
Large	31.2	960
Extra Large	31.3	970
ASTM 3578:05	18 MPa minimum	650% minimum

AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	28.0	900
Small	28.0	920
Medium	27.5	910
Large	27.3	890
Extra Large	27.9	900
ASTM 3578:05	14 MPa minimum	500% minimum

DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
6.5 Aurelia Dimensions ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
7 Aurelia Dimensions ASTM 3578:05	9.4 / 240 8.7 / 220 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
7.5 Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
8 Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum
8.5 Aurelia Dimensions ASTM 3578:05	9.4 / 240 9.1 / 230 minimum	5.0 / 0.12 3.2 / 0.08 minimum	6.0 / 0.15 3.2 / 0.08 minimum

ITEM NUMBER

6.5	7	7.5	8	8.5
69925	69926	69927	69928	69929

Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.



Authorized Aurelia Distributor:

CAUTION STATEMENT: This product contains natural rubber latex which may cause allergic reactions. Safe use of this glove by or on latex sensitized individuals has not been established.

AURELIA Div. Supermax Inc.

2225 White Oak Circle, Aurora, IL 60502 | Toll Free: 1-877-AURELIA | Tel: 1-630-898-8886 | Fax: 1-630-898-8855
Web: www.aureliagloves.com | Email: info@aureliagloves.com