

# Dare to Compare.

## Aurelia® Quest™ White Nitrile Exam Gloves

- Latex Comfort
- Latex Fit
- Latex Cost
- Nitrile Durability
- Nitrile-Latex Free

## Latex Gloves

- Latex Comfort
- Latex Fit
- Latex Cost
- Nitrile Durability
- Nitrile-Latex Free


## Aurelia® Quest™

### The Finest Nitrile Glove Ever Made.

At last, a nitrile exam glove that fits and feels like Latex, with all the other Nitrile benefits like superior durability, and Non-Latex allergenic. Aurelia Quest delivers the tactile sensitivity, comfort and fit of Latex. Even in the fingertips! Ultra-sheer yet durable. The most comfortable Nitrile ever, in a clinical white color! Try for yourself! Call or email:

[QuestSampleRequest@AureliaGloves.com](mailto:QuestSampleRequest@AureliaGloves.com)

1-877-AURELIA (1-877-287-3542) [www.aureliagloves.com](http://www.aureliagloves.com)



**AURELIA®**  
Sign of Comfort™

# Quest

## Product Description

Aurelia® Quest™ White Nitrile Powder Free Examination Gloves are made from 100% Nitrile. This New Generation Nitrile gloves are super soft & lightweight which contour to your hands for superior comfort. Each batch of gloves are air tested to detect pinholes. Aurelia gloves conform and are tested to the highest international standards. For more information, visit [www.aureliagloves.com](http://www.aureliagloves.com)



- 100% Latex Free • White Color
- New Generation of Nitrile Glove for superior comfort
- Finger Textured for precise gripping
- Greater chemical resistance compared to other materials
- Non Sterile, Ambidextrous, Single Use
- 100 gloves per box

## Quest® Dimensions™ Specifications

### TENSILE TEST RESULTS (unaged)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	680
Small	31.1	680
Medium	31.0	680
Large	31.2	680
Extra Large	31.3	680
ASTM D6319	14 MPa minimum	500% minimum

### AGED TENSILE TEST RESULTS (7 Days at 70° C)

	Tensile Strength (MPa)	Elongation at break (%)
Extra Small	30.0	620
Small	30.0	620
Medium	30.0	620
Large	30.0	620
Extra Large	30.0	620
ASTM D6319	14 MPa minimum	620% minimum

### DIMENSIONS

	Length (inches/mm)	Thickness (Single-wall)	
		Palm (mils/mm)	Finger (mils/mm)
<b>Extra Small</b> Aurelia Quest ASTM D6319	9.4 / 240 8.7 / 220 minimum	3.0 / 0.08 2.0 / 0.05 minimum	4.0 / 0.10 2.0 / 0.05 minimum
<b>Small</b> Aurelia Quest ASTM D6319	9.4 / 240 8.7 / 220 minimum	3.0 / 0.08 2.0 / 0.05 minimum	4.0 / 0.10 2.0 / 0.05 minimum
<b>Medium</b> Aurelia Quest ASTM D6319	9.4 / 240 8.7 / 220 minimum	3.0 / 0.08 2.0 / 0.05 minimum	4.0 / 0.10 2.0 / 0.05 minimum
<b>Large</b> Aurelia Quest ASTM D6319	9.4 / 240 8.7 / 220 minimum	3.0 / 0.08 2.0 / 0.05 minimum	4.0 / 0.10 2.0 / 0.05 minimum
<b>Extra Large</b> Aurelia Quest ASTM D6319	9.4 / 240 8.7 / 220 minimum	3.0 / 0.08 2.0 / 0.05 minimum	4.0 / 0.10 2.0 / 0.05 minimum


## Ordering Information

To order from your preferred distributor, specify size and item number. Packing: 100 pieces per dispenser box, 10 dispenser boxes per shipping carton.

### ITEM NUMBER

Extra Small	Small	Medium	Large	Extra Large
92895	92896	92897	92898	92899

Authorized Aurelia Distributor:



**AURELIA**  
Sign of Comfort™